

Ephys pipeline overview

Folder structure

- Y:\Data\TDTtanks\Monkey_phys contains the Raw TDT data
- Y:\Data\Sortcodes\Monkey_phys contains sorting related files
- Y:\Data\Monkey contains the behavioral data
- Y:\Data\Monkey_phys_mat_from_TDT contains recorded data, arranged in the same trial structure as the behavioral data
- Y:\Data\Monkey_phys_combined_monkeypsych_TDT contains trial structures with combined behavioral and ephys data
- Y:\Data\All_phys_preprocessing_log\Monkey_phys contains automatically created log files of all preprocessing performed
- \Dropbox\DAG\phys\Monkey_phys_dpz contains
 1. Electrode_depths matfile
 2. Same_cells matfile
 3. plx files excel table
 4. Sorted_neurons excel table

General Workflow

- Update Electrode depths matfile. Only block/channel combinations defined here will be processed (?)
- (Optional - Synapse only) Reformat Synapse tank folder format to match previous conventions
- (Optional - WC only) Run waveclus automatic sorting (Channels with electrodes in the same depth across blocks will be combined across blocks)
- (Optional - WC only) manually refine WC sorting
- (Optional - Plexon only) Create Plexon sortcodes
- (Optional - Plexon only) Manual Plexon sorting
- Update plx files excel table
- Combine TDT and behavioral data
 1. Create TDT trial structure
 2. Combine
- Update Same_cells matfile
- Update sorted_neurons excel table

Most of the steps above are automatized, check out phys_gui page for details.

Synchronization

- TDT data is stored as a continuous data (starting from start of recording)
- behavioral data is stored per trial (starting from first trial)

For combining both, we first convert TDT data to trial format (function TDT_trial_struct.m) preferentially using epochs store (Tnum and SVal) information.

Note: There is a flag for using the continuous state information stream (stream_state_info), but it is permanently turned off, as the epoch information is already extracted from the continuous state information within the TDT circuit.

Importantly, trial start reference is STATE 2 (fixation acquisition). Trial end is STATE 2 of the next trial. Therefore, ITI between trial 1 and 2 as well as STATE 1 (initiation) of trial 2 are found at the end of trial 1.

This also means that data recorded before 1st trial's state 2 can not be stored in the trial structure and is saved in a separate variable (First_trial_INI) instead.

Associated code

From:

<http://dag.dokuwiki.dpz.lokal/> - **DAG wiki**

Permanent link:

http://dag.dokuwiki.dpz.lokal/doku.php?id=ephys_pipeline:1_pipeline_overview&rev=1639495076

Last update: **2022/12/29 07:15**

