

# Spike analysis pipeline overview

## Folder structure

Inputs	Info
\Dropbox\DAG\phys\Monkey_phys_dpz\Mon_sorted_neurons.xls	electrode locations, unit/site IDs,SNR/Stability/single ratings
\GitHub\Settings\spike_analysis\project	project and version specific settings (keys)
Y:\Data\Monkey_phys_combined_monkeypsych_TDT	raw data

Outputs in Y:\Projects\project\ephys\version:

File	Info
population_monkey_session.mat	spike data sorted by unit for each session
sites_monkey_session.mat	LFP data sorted by site for each session
by_block_monkey_session.mat	body signals sorted by block for each session
Mon_sorted_neurons.xls	copies of the used sorted neurons tables from dropbox
keys_Monkey.mat	copy of the keys used when running ph_initiation
tuning_table_combined.mat	Anova results stored in a table
tuning_table_combined_CI.mat	Anova results restructured according to contra/ipsi definitions
tuning_table_combined.xls	simplified excel table
seed.mat	a saved seed to make randomizations reproducible
Subfolder	Info
spike_shapes	Spike shapes, firing rates over time, and ISI plots
single_cell_examples	Single cell plots
cell_counts	Anova results as pie plots
scatter	Anova results as scatter plots (one column versus another)
population_analysis	population PSTHs
response timing	tuning over time plots
...	

## General Workflow

- ph\_initiation(project,{version1,version2,...})
  1. read in general settings → project settings → version settings
  2. loop per monkey: ph\_session\_processing (core function)
    - read in data from Y:\Data\Monkey\_phys\_combined\_monkeypsych\_TDT
    - run monkeypsych\_analyze for saccade detection etc.
    - run ph\_run\_state\_alignment\_per\_trial
      - takes over relevant trial and state (event) information from monkeypsych\_analyze output (TDT\_states and saccade/reach initiation)
      - combines it with Sorted Neuron table information
      - !! Copies last 1 second of spikes to beginning of next trial
      - !! Cuts and appends last 1 second of streams to beginning of next trial
      - !! excludes trials without physiology data And/or NOT matching condition

- resort data by unit/site/block
  - plot waveforms/ISI/FR\_across time per unit
  - exclude units dependent on excel entries (SNR/stability/single rating)
  - plot waveforms/ISI/FR\_across time per unit again for remaining units
  - run ANOVAs (ph\_ANOVAS) and create tuning table
  - create single cell plots (ph\_plot\_unit\_per\_condition)
  - save files per session
3. format tuning table (create tuning\_table\_combined\_CI.mat and tuning\_table\_combined.xls)
  4. ph\_get\_filelist for crossreferencing with behavior
  5. ph\_initiate\_population\_analysis also runs population analysis

## Synchronization

- Very important: In additi

## Associated code

hmm

From:

<http://dag.dokuwiki.dpz.lokal/> - DAG wiki

Permanent link:

[http://dag.dokuwiki.dpz.lokal/doku.php?id=spike\\_analysis\\_pipeline:1\\_pipeline\\_overview&rev=1672298151](http://dag.dokuwiki.dpz.lokal/doku.php?id=spike_analysis_pipeline:1_pipeline_overview&rev=1672298151)

Last update: 2022/12/29 07:15

