

# Spike analysis pipeline overview

## Folder structur

Inputs	Info
\\Dropbox\\DAG\\phys\\Monkey_phys_dpz\\Mon_sorted_neurons.xls	electrode locations, unit/site IDs, SNR/Stability/single ratings
\\GitHub\\Settings\\spike_analysis\\project	project and version specific settings (keys)
Y:\\Data\\Monkey_phys_combined_monkeypsych_TDT	raw data

Outputs in Y:\\Projects\\project\\ephys\\version:

File	Info
population_monkey_session.mat	spike data sorted by unit for each session
sites_monkey_session.mat	LFP data sorted by site for each session
by_block_monkey_session.mat	body signals sorted by block for each session
Mon_sorted_neurons.xls	copies of the used sorted neurons tables from dropbox
keys_Monkey.mat	copy of the keys used when running ph_initiation
tuning_table_combined.mat	Anova results stored in a table
tuning_table_combined_CI.mat	Anova results restructured according to contra/ipsi definitions
tuning_table_combined.xls	simplified excel table
seed.mat	a saved seed to make randomizations reproducible
Subfolder	Info
spike_shapes	Spike shapes, firing rates over time, and ISI plots
single_cell_examples	Single cell plots
cell_counts	Anova results as pie plots
scatter	Anova results as scatter plots (one column versus another)
population_analysis	population PSTHs
response timing	tuning over time plots
...	

## General Workflow

- `ph_initiation(project, {version1, version2, ...})`
  1. read in general settings → project settings → version settings
  2. loop per monkey: `ph_session_processing` (core function)
    - read in data from Y:\\Data\\Monkey\_phys\_combined\_monkeypsych\_TDT
    - run `monkeypsych_analyze` for saccade detection etc.
    - run `ph_run_state_alignment_per_trial`
      - takes over relevant trial and state (event) information from `monkeypsych_analyze` output (TDT\_states and saccade/reach initiation)
      - combines it with Sorted Neuron table information
      - !! Copies last 1 second of spikes to beginning of next trial
      - !! Cuts and appends last 1 second of streams to beginning of next trial
      - !! excludes trials without physiology data And/or NOT matching condition

- resort data by unit/site/block
  - plot waveforms/ISI/FR\_across time per unit
  - exclude units dependent on excel entries (SNR/stability/single rating)
  - plot waveforms/ISI/FR\_across time per unit again for remaining units
  - run ANOVAs (ph\_ANOVAS) and create tuning table
  - create single cell plots (ph\_plot\_unit\_per\_condition)
  - save files per session
3. format tuning table (create tuning\_table\_combined\_CI.mat and tuning\_table\_combined.xls)
  4. ph\_get\_filelist for crossreferencing with behavior
  5. ph\_initiate\_population\_analysis also runs population analysis

## Synchronization

- Very important: ... INDEED :)

## Associated code

...

From:  
<http://dag.dokuwiki.dpz.lokal/> - DAG wiki

Permanent link:  
[http://dag.dokuwiki.dpz.lokal/doku.php?id=spike\\_analysis\\_pipeline:1\\_pipeline\\_overview&rev=1699894423](http://dag.dokuwiki.dpz.lokal/doku.php?id=spike_analysis_pipeline:1_pipeline_overview&rev=1699894423)

Last update: 2023/11/13 16:53

