

Spike analysis pipeline overview

Folder structure

- Y:\Projects\project\ephys\version main folder contains:
 1. population_monkey_session.mat : spike data sorted by unit for each session
 2. sites_monkey_session.mat : LFP data sorted by site for each session
 3. by_block_monkey_session.mat : body signals sorted by block for each session
 4. Mon_sorted_neurons.xls : copies of the used sorted neurons tables from dropbox
 5. keys_Monkey.mat : copy of the keys used when running ph_initiation
 6. tuning_table_combined.mat : Anova results stored in a table
 7. tuning_table_combined_CI.mat : Anova results restructured according to contra/ipsi definitions
 8. tuning_table_combined.xls : simplified excel table
 9. seed.mat : a saved seed to make randomizations reproducible
- Subfolders:
 1. spike_shapes : Spike shapes, firing rates over time, and ISI plots
 2. single_cell_examples : Single cell plots
 3. cell_counts : Anova results as pie plots
 4. scatter : Anova results as scatter plots (one column versus another)
 5. population_analysis : population PSTHs
 6. response timing : tuning over time plots
 7. ...
- \Dropbox\DAG\phys\Monkey_phys_dpz\Sorted_neurons excel table

General Workflow

- ph_initiation(project, {version1, version2, ...})
 1. loop per monkey
 2. read in general settings → project settings → version settings
 3. run ph_session_processing (core function)
 - read in data from Y:\Data\Monkey_phys_combined_monkeypsych_TDT
 - run monkeypsych_analyze for saccade detection etc.
 - run ph_run_state_alignment_per_trial
 1. takes over relevant trial and state (event) information
 2. combines it with Sorted Neuron table information
 3. !! Copies last 1 second of spikes to beginning of next trial
 4. !! Cuts and appends last 1 second of streams to beginning of next trial
 5. !! excludes trials without physiology data And/or NOT matching condition
 - resort data by unit/site/block
 - plot waveforms/ISI/FR_across time per unit
 - exclude units dependent on excel entries (SNR/stability/single rating)
 - plot waveforms/ISI/FR_across time per unit again for remaining units
 - run ANOVAs (ph_ANOVAS)
 - create single cell plots (ph_plot_unit_per_condition)
 - save files per session
 4. format tuning table (create tuning_table_combined_CI.mat and

tuning_table_combined.xls)

5. ph_get_filelist for crossreferencing with behavior

6. ph_initiate_population_analysis also runs population analysis

Synchronization

- Very important: In additi

Associated code

hmm

From:
<http://dag.dokuwiki.dpz.lokal/> - **DAG wiki**

Permanent link:
http://dag.dokuwiki.dpz.lokal/doku.php?id=spike_analysis_pipeline:pipeline_overview&rev=1641838055

Last update: **2022/12/29 07:15**

