

Spike analysis pipeline overview

Folder structure

Inputs	Info
\\Dropbox\\DAG\\phys\\Monkey_phys_dpz\\Mon_sorted_neurons.xls	electrode locations, unit/site IDs, SNR/Stability/single ratings
\\GitHub\\Settings\\spike_analysis\\project	project and version specific settings (keys)
Y:\\Data\\Monkey_phys_combined_monkeypsych_TDT	raw data

Outputs in Y:\\Projects\\project\\ephys\\version:

File	Info
population_monkey_session.mat	spike data sorted by unit for each session
sites_monkey_session.mat	LFP data sorted by site for each session
by_block_monkey_session.mat	body signals sorted by block for each session
Mon_sorted_neurons.xls	copies of the used sorted neurons tables from dropbox
keys_Monkey.mat	copy of the keys used when running ph_initiation
tuning_table_combined.mat	Anova results stored in a table
tuning_table_combined_CI.mat	Anova results restructured according to contra/ipsi definitions
tuning_table_combined.xls	simplified excel table
seed.mat	a saved seed to make randomizations reproducible
Subfolder	Info
spike_shapes	Spike shapes, firing rates over time, and ISI plots
single_cell_examples	Single cell plots
cell_counts	Anova results as pie plots
scatter	Anova results as scatter plots (one column versus another)
population_analysis	population PSTHs
response timing	tuning over time plots
...	

General Workflow

- `ph_initiation(project, {version1, version2, ...})`
 1. read in general settings → project settings → version settings
 2. loop per monkey: `ph_session_processing` (core function)
 - read in data from Y:\\Data\\Monkey_phys_combined_monkeypsych_TDT
 - run `monkeypsych_analyze` for saccade detection etc.
 - run `ph_run_state_alignment_per_trial`
 - takes over relevant trial and state (event) information from `monkeypsych_analyze` output (TDT_states and saccade/reach initiation)
 - combines it with Sorted Neuron table information
 - !! Copies last 1 second of spikes to beginning of next trial
 - !! Cuts and appends last 1 second of streams to beginning of next trial
 - !! excludes trials without physiology data And/or NOT matching condition

- resort data by unit/site/block
 - plot waveforms/ISI/FR_across time per unit
 - exclude units dependent on excel entries (SNR/stability/single rating)
 - plot waveforms/ISI/FR_across time per unit again for remaining units
 - run ANOVAs (ph_ANOVAS) and create tuning table
 - create single cell plots (ph_plot_unit_per_condition)
 - save files per session
3. format tuning table (create tuning_table_combined_CI.mat and tuning_table_combined.xls)
 4. ph_get_filelist for crossreferencing with behavior
 5. ph_initiate_population_analysis also runs population analysis

Synchronization

- Very important: In additi

Associated code

hmm

From:
<http://dag.dokuwiki.dpz.lokal/> - **DAG wiki**

Permanent link:
http://dag.dokuwiki.dpz.lokal/doku.php?id=spike_analysis_pipeline:pipeline_overview&rev=1652798354

Last update: **2022/12/29 07:15**

