

# AFNI analysis

## Converting individual (ACPC) monkey brains to NMT

```
3dAllineate -lDmatrix_apply curius_anat_composite_linear_to_template.1D -
prefix CU_NMT -source CUclean_20130627_ACPC.nii

3dNwarpApply -nwarp curius_anat_shft_WARP.nii.gz -source CU_NMT+tlrc.BRIK

3dAFNItoNIFTI CU_NMT_Nwarp+tlrc.BRIK
```

## Converting VMPs from original ACPC space to NMT

E.g. Y:\MRI\Curius\2NMT\vmp

Original vmp: dPul\_p.vmp, “nominal” res 2 mm, actual 1 mm

The approach works with multiple maps within vmp

```
vmp.ExportNifti('dPul_p.nii'); % or NE GUI
ne_change_nii_header_resolution('dPul_p.nii',1);
```

```
igor@alex-Precision-T3610:~/VMP_warps$
cp /mnt/KognitiveNeurowissenschaften/DAG/MRI/Curius/2NMT/vmp/dPul_p.nii .
3dAllineate -master CU_NMT_Nwarp.nii -newgrid 1 -lDmatrix_apply
curius_anat_composite_linear_to_template.1D -prefix dPul_p -source
dPul_p.nii
3dNwarpApply -nwarp curius_anat_shft_WARP.nii.gz -source dPul_p+tlrc.BRIK
3dAFNItoNIFTI dPul_p_Nwarp+tlrc.BRIK
sudo cp dPul_p_Nwarp.nii
/mnt/KognitiveNeurowissenschaften/DAG/MRI/Curius/2NMT/vmp/
```

```
ne_change_nii_header_resolution('dPul_p_Nwarp.nii',2);
ne_afni_nii2vmp('dPul_p_Nwarp.nii','t',2,'','dPul_p.vmp');
```

## Converting electrode localization prob. maps VMPs from original ACPC space to NMT

E.g. Y:\MRI\Curius\2NMT\electrode\_localization

Original vmp: el\_dPul\_p.vmp, “nominal” res 1 mm, actual 0.5 mm

```
ne_vmp2vmr2nifti('el_dPul_p.vmp','CUclean_20130627_ACPC_1mm.vmr');
```

```
ne_change_nii_header_resolution('el_dPul_p.nii',0.5);
```

```
igor@alex-Precision-T3610:~/VMP_warps$  
cp /mnt/KognitiveNeurowissenschaften/DAG/MRI/Curius/2NMT/vmp/el_dPul_p.nii .  
3dAllineate -master CU_NMT_Nwarp.nii -ldmatrix_apply  
curius_anat_composite_linear_to_template.1D -prefix el_dPul_p -source  
el_dPul_p.nii  
3dNwarpApply -nwarp curius_anat_shft_WARP.nii.gz -source el_dPul_p+tlrc.BRIK  
3dAFNItoNIFTI el_dPul_p_Nwarp+tlrc.BRIK  
sudo cp el_dPul_p_Nwarp.nii  
/mnt/KognitiveNeurowissenschaften/DAG/MRI/Curius/2NMT/vmp/
```

```
ne_change_nii_header_resolution('el_dPul_p_Nwarp.nii',1);  
ne_afni_nii2vmp('el_dPul_p_Nwarp.nii','t',1);
```

From:

<http://dag.dokuwiki.dpz.lokal/> - DAG wiki

Permanent link:

<http://dag.dokuwiki.dpz.lokal/doku.php?id=analysis:fmri:afni&rev=1615292566>

Last update: 2022/12/29 07:15

